RibLoc® U Plus Chest Wall Plating System

| | | |
|----|-------------------------------|-----|
| US | Instructions for use | 2 |
| DK | Brugsanvisning..... | 12 |
| DE | Gebrauchsanweisung | 23 |
| EL | Οδηγίες χρήσης..... | 35 |
| EN | Instructions for use | 48 |
| ES | Instrucciones de uso | 57 |
| FI | Käyttöohjeet | 69 |
| FR | Mode d'emploi..... | 80 |
| ID | Istruzioni per l'uso | 92 |
| IT | Bruksanvisning..... | 104 |
| NO | Gebrauksinstructies..... | 115 |
| NL | Instruções de utilização..... | 127 |
| PT | Bruksanvisning..... | 139 |
| SV | Kullanım talimatları..... | 150 |

Instructions for use

RibLoc® U Plus Chest Plating System

These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals. The US instructions are intended for users in the United States and its territories.

Rx only

DESCRIPTION

The RibLoc U Plus System of bone plates, screws and instruments are designed to provide fixation for fractures, fusions, and osteotomies of the ribs and reconstruction of the chest wall.

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CONTRAINDICATIONS

- Active or latent infection.
- Sepsis.
- Insufficient quantity or quality of bone, osteoporosis.
- Patients with confirmed material sensitivity.
- Patients who are unwilling or incapable of following post-operative care instructions.
- These devices are not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

WARNINGS & PRECAUTIONS

Warning:

- For safe and effective use of this implant, the surgeon must be thoroughly familiar with the implant, the methods of application, the instruments and the recommended surgical technique for this device.
- Surgeons must carefully consider the likelihood of healing being achieved when plating fractures, osteotomies, or reconstructions of the chest wall. This system is only designed to withstand loading during a reasonable healing time period and is not intended to be a permanent prosthesis.
- Surgeons must consider a possible need for emergent reentry, such as sternotomy, before plating the sternum.
- Improper insertion of the device during implantation can increase the possibility of loosening or migration.
- A poorly contoured plate may result in an abnormal load on bones or may result in patient discomfort.
- Device damage or breakage can occur when the implant is subjected to increased loading associated with trauma, delayed union, nonunion or incomplete healing. Device breakage could lead to additional surgery and device removal.
- The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or continuous load bearing past

the average healing time (6-8 weeks), particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete bone healing.

- The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail.
- As with any surgical implantation there is a possibility of nerve, bone or soft tissue damage related to either surgical trauma or the presence of the implant.

Caution:

- An implant shall never be reused. Previous stresses may have created imperfections which can lead to a device failure.
- Extreme or repeated bending of the implants can cause stresses that may lead to premature device failure.
- Ensure that RibLoc plates and U-clips are not damaged prior to installation and use care when handling to prevent U-clip deformation.
- Over-compression of the U-Clip during implant installation can damage the bone.
- Instruments should be inspected for wear or damage prior to usage.
- The drill bit shall be discarded after each surgery since after normal use the drill bit can become too dull to perform as intended.
- If a plate cutter is used on the implant take necessary precautions, a sharp edge may have been created.
- During use of a drill, cutting or installing plates and inserting screws, take necessary precautions when in close proximity to sharp edges and points, and be aware that metal debris/fragments can be generated. Remove any observed debris/fragments from the surgical field with suction or manually and dispose of appropriately.
- Protect implants against scratching and nicking prior to and during use. Prevent hexalobe drivers, drill bits and instruments to come into inadvertent contact with plates as such stress concentrations can lead to device failure.
- The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

ADVERSE EFFECTS**Possible adverse effects include:**

- Pain, discomfort, or abnormal sensations.
- Nerve or soft tissue damage.
- Necrosis of bone or tissue.
- Bone resorption.
- Inadequate healing from the presence of an implant or due to surgical trauma
- Implant fracture, migration and/or loosening may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion. Implant migration and/or loosening may occur.
- Metal sensitivity, histological, allergic or adverse foreign body reaction resulting from implantation of a foreign material. Consult our document "Metal Sensitivity Statement" at www.acumed.net/ifu.

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The Surgical Technique also includes information about the installation and calibration of the device.

Important: Surgical Techniques may contain important safety information.

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MRI SAFETY INFORMATION

Many Acumed implants have been evaluated for safety in the MR environment and are MR conditional, with exception for sternum plate which is MRI Unsafe. Consult our publication "Acumed Implants in the MR Environment" at www.acumed.net/ifu for more information.

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IMPLANTS

MATERIALS

- The implant is made of commercially pure Titanium per ASTM F67.
- Consult our document "Metal Sensitivity Statement" at www.acumed.net/ifu for the chemical composition of Acumed metal implants.

- Chemical characterization data confirmed that all identified chemicals and elements were present at levels below 0.1 w/w%. Toxicological assessment further showed that all identified values were below tolerable exposure levels and are biologically safe.

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MULTIPLE USE and SINGLE USE

- Instruments are intended for multiple use unless identified on the label for single use only.
- Single use instruments are intended to be disposed after use on a single patient during a single procedure.
- Do not reuse single use instruments as this may increase the risks of failure and cross-contamination.

- Multiple use instruments are only intended for use on a single patient and a single procedure before requiring processing.
- Multiple (limited) use instruments have a limited lifespan. Immediately replace any multiple use instrument if performance becomes inadequate.
- Multiple use instruments potentially contaminated with transmissible spongiform encephalopathy (TSE) agents shall not be processed for reuse.

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IMPORTANT

- Promptly perform the processing steps to limit microbial growth and maximize the effectiveness of sterilization.
- Inspect implants for contamination by blood or tissue and dispose of them when found. Do not process contaminated implants.
- Prevent instrument corrosion by minimizing contact with solutions containing iodine, chlorine, and saline or other metal salts.
- Prevent damage to the protective anodization layer on aluminum instruments by avoiding contact with solutions < 4 pH and > 9 pH, especially if they contain sodium carbonate or sodium hydroxide.
- Repeated processing of anodized metals may cause colors to fade but this does not affect the function of the device.
- Avoid cleaning agents containing aldehydes since they can denature and coagulate proteins (fixation).
- Enzymatic detergents are well suited for loosening protein-based contamination.
 - Use a neutral pH enzymatic detergent.
 - Use a low foaming solution to allow visibility of the device during cleaning.
- Closely follow the manufacturer's instructions for the safety, storage, mixing, water quality, exposure time, temperature, replacement, and disposal of cleaning agents.
- Devices potentially contaminated with transmissible spongiform encephalopathy (TSE) agents shall not be processed or reused. These processing instructions are not suitable for inactivation of TSE agents. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling and disposal of devices potentially contaminated with TSE agents.
- Utility water: Refer to AAMI TIR34* when instructed to use utility water. Utility water is typically municipal or tap water but may require additional treatment to be suitable for use.

- Critical water: Refer to AAMI TIR34* when instructed to use critical water. Critical water is highly treated and has very low organic and inorganic content with an endotoxin level under 10
- EU/mL. Suitable water may also be specified in national pharmacopeias, national standards, and hospital protocols.

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MANUAL CLEANING

1. Rinse the contaminated instruments under running cold utility water to reduce heavy surface contamination.
2. Dispose of any used instruments intended for single use only.
3. Place the contaminated instruments in enzymatic solution* until completely submerged to minimize the spraying of solution.
4. Actuate all moveable parts to allow detergent to contact all surfaces
5. Soak for a minimum of ten (10) minutes.
6. Scrub the instruments using a soft-bristled brush to remove all visible debris. Do not use stainless steel or other abrasives as these may damage the surface.
 - When possible, scrub the instruments when totally submerged to minimize the spraying of fluid.
7. Some instruments may require special consideration:
 - Clean the instruments with all parts loosened. Clean the instruments disassembled if they are designed to be taken apart.
 - Use a water jet to flood cleaning solution into challenging areas, such as mating surfaces, springs, coils, cannulations, blind holes, flutes, cutting teeth, and flexible parts to flush out any trapped soil.
 - Operate movable parts and rotate (as necessary) while scrubbing to ensure that all crevices are accessible.
 - Carefully clean cannulated parts and challenging areas using an appropriately sized brush.
 - Optionally sonicate for 10 to 15 minutes using a fresh, neutral pH ultrasonic cleaning solution. Follow the ultrasonic cleaner and detergent manufacturer's instructions. Important! Any previous surface damage may increase due to ultrasonic cleaning.
8. Perform an initial rinse for at least 3 minutes using clean, soft, utility water in the temperature range of 25°C to 35°C (77°F to 95°F) to remove all signs of contamination and cleaning agent.
 - Actuate all movable parts.
 - Flush out cannulations and complex mechanisms.
9. Repeat the previous processing steps if visible residue remains present.
10. Perform a final rinse for at least 1 minute using critical water to displace minerals and other impurities found in utility water. Do not use saline solutions for final rinsing because they may interfere with disinfection and sterilization.
 - Actuate all moving parts
 - Pay particular attention to cannulations and blind holes as well as hinges and joints between mating surfaces.
Rinse cannulations at least three times with a syringe (volume 1-50ml).
11. Remove excess moisture from the instruments using a clean, absorbent, non-shedding wipe. 12. Allow the instruments to thoroughly dry. Any moisture may affect sterilization and devices may remain wet after the drying period.

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STERILIZATION

- Perform sterilization using a dynamic-air-removal (prevacuum) autoclave. Gravity displacement sterilization is not recommended. Immediate use (flash) sterilization is not recommended.
 - Ensure the sterilizer’s maximum load limit is not exceeded when sterilizing multiple sets or devices.
 - Do not stack containers as this might prevent the penetration of steam and inhibit drying.
 - Refer to the sterilizer manufacturer's instructions and ensure proper installation, calibration, use, and ongoing maintenance.
 - The sterilized items should be allowed to cool to room temperature before handling. This allows for safe handling and preventing condensation.
 - Follow current industry best practice guidelines such as ANSI/AAMI ST79:2017*.
- * Association for the Advancement of Medical Instrumentation (AAMI). Comprehensive guide to steam sterilization and sterility assurance in health care facilities. AAMI ST79:2017. Arlington, VA.
- The following table shows the minimum parameters validated* to achieve a required Sterility Assurance Level (SAL) of 10⁻⁶ for the system.

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STORAGE CONDITIONS

STORAGE OF PACKAGED PRODUCT

- Items should be stored under controlled conditions in a manner that minimizes the potential for contamination per ANSI/AAMI ST79:2017.

- Refer to sterilization wrap or rigid container manufacturer's IFU for limits on sterile product storage time and storage requirements for temperature and humidity.

STORAGE OF STERILE PRODUCT

- Final packaged product should be stored at room temperature and protected from direct sunlight, pests, and high humidity.

Symbols Glossary

| Symbol | Description | ISO 15223-1 |
|---|--|---------------------|
|  | Consult the electronic instructions for use (eIFU) at www.acumed.net/ifu | 5.4.3 |
|  | Caution | 5.4.4 |
|  | Non-sterile | 5.2.7 |
|  | Catalogue number | 5.1.6 |
|  | Batch code | 5.1.5 |
|  | Authorized representative in the European Community / European Union | 5.1.2 |
|  | Medical device | 5.7.7 |
|  | Manufacturer | 5.1.1 |
|  | Date of manufacture | 5.1.3 |
|  | Do not re-use | 5.4.2 |
|  | Do not use if package is damaged and consult instructions for use / do not use if the product sterile barrier system or its packaging is compromised | 5.2.8 |
| Rx Only | Caution: U.S. federal law restricts this device to sale by or on the order of a physician. | U.S. 21 CFR 801.109 |
|  | The reticle is a registered trademark of Acumed. It may appear alone or with the Acumed name. | |
|  | CE marking of conformity, Article 17 of EU Directive 93/42/EEC or Article 20 of Regulation (EU) 2017/745. CE marking may be accompanied by the identification number of the notified body responsible for conformity assessment. | |
|  | MR Conditional: An item with demonstrated safety in the MR environment within defined conditions. | |

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RibLoc U Plus Chest Plating System

EN

These instructions are intended for the Operating Surgeon and supporting Healthcare Professionals. The EN instructions are intended for users in English speaking countries outside the United States and its territories.

DESCRIPTION

The RibLoc U Plus System of bone plates, screws and accessories are designed to provide fixation for fractures, fusions, and osteotomies of the ribs and reconstruction of the chest wall.

INDICATIONS FOR USE

The RibLoc U Plus System is intended to stabilize and provide fixation for fractures, fusions, and osteotomies of the ribs and for reconstructions of the chest wall, and sternum.

CONTRAINDICATIONS

- Active or latent infection.
- Sepsis.
- Insufficient quantity or quality of bone/ soft tissue and material sensitivity If material sensitivity is suspected, tests should be performed prior to implantation.
- Patients who are unwilling or incapable of following post-operative care instructions.
- These devices are not intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

WARNINGS / RESIDUAL RISKS & PRECAUTIONS

Warnings & residual Risks include:

- For safe and effective use of this implant, the surgeon must be thoroughly familiar with the implant, the methods of application, the instruments and the recommended surgical technique for this device.
- Surgeons must carefully consider the likelihood of healing being achieved when plating fractures, osteotomies, or reconstructions of the chest wall. This system is only designed to withstand loading during a reasonable healing time period and is not intended to be a permanent prosthesis.
- Surgeons must consider a possible need for emergent reentry, such as sternotomy, before plating the sternum.
- Improper insertion of the device during implantation can increase the possibility of loosening or migration.
- A poorly contoured plate may result in an abnormal load on bones or may result in patient discomfort.
- Device damage or breakage can occur when the implant is subjected to increased loading associated with trauma, delayed union, nonunion or incomplete healing. Device breakage could lead to additional surgery and device removal.
- The patient must be cautioned, preferably in writing, about the use, limitations, and possible adverse effects of this implant. These cautions include the possibility of the device or treatment failing as a result of loose fixation and/or loosening, stress, excessive activity, or continuous load bearing past the

average healing time (6-8 weeks), particularly if the implant experiences increased loads due to delayed union, nonunion, or incomplete bone healing.

- The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail.
- As with any surgical implantation there is a possibility of nerve, bone or soft tissue damage related to either surgical trauma or the presence of the implant.
- The implant system has not been evaluated for use in pectus deform.

Precautions:

- An implant shall never be reused. Previous stresses may have created imperfections which can lead to a device failure.
- Extreme or repeated bending of the implants can cause stresses that may lead to premature device failure.
- Ensure that RibLoc plates and U-clips are not damaged prior to installation and use care when handling to prevent U-clip deformation.
- Over-compression of the U-Clip during implant installation can damage the bone.
- Instruments should be inspected for wear or damage prior to usage.
- The drill bit shall be discarded after each surgery since after normal use the drill bit can become too dull to perform as intended.
- If a plate cutter is used on the implant take necessary precautions, a sharp edge may have been created.
- During use of a drill, cutting or installing plates and inserting screws, take necessary precautions when in close proximity to sharp edges and points, and be aware that metal debris/fragments can be generated. Remove any observed debris/fragments from the surgical field with suction or manually and dispose of approximately.
- Protect implants against scratching and nicking prior to and during use. Prevent hexalobe drivers, drill bits and instruments to come into inadvertent contact with plates as such stress concentrations can lead to device failure.
- The benefits from implant surgery may not meet the patient's expectations or may deteriorate with time necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are not uncommon.

ADVERSE EFFECTS

Possible adverse effects include:

- Pain, discomfort, or abnormal sensations due to the presence of an implant.
- Nerve or soft tissue damage, necrosis of bone or bone resorption, necrosis of the tissue or inadequate healing from the presence of an implant or due to surgical trauma.
- Implant fracture, migration and/or loosening may occur due to excessive activity, prolonged loading upon the device, incomplete healing, or excessive force exerted on the implant during insertion.
- Metal sensitivity, histological, allergic or adverse foreign body reaction resulting from implantation of a foreign material. Consult our document "Metal Sensitivity Statement" at www.acumed.net/ifu.

Important: Any serious incident that has occurred in relation to the device should be reported to Acumed through complaints@acumed.net.

TARGET POPULATION

The implant is used to treat adult patient populations who suffer from fractures, fusions, and osteotomies of the ribs and for reconstructions of the chest wall, and sternum.

SURGICAL TECHNIQUE

Acumed offers one or more Surgical Techniques to promote the safe and effective use of this system. Consult our Surgical Techniques at www.acumed.net.

Important: Surgical techniques may contain important safety information.

Important: The instruments and implants in this system are intended to be used by suitably trained and qualified surgeons in a hospital operating room setting. Before treatment, the surgeon is advised to read and fully understand all instructions and communicate to the patient any relevant medical information provided therein, including the use, limitations, risks (safety communications), and possible adverse effects of the proposed treatment.

Consult the most recent versions of the Instructions for Use and Surgical Techniques as they are subject to change. Contact Acumed or an authorized agent to request any additional information.

SYSTEM COMPATIBILITY

- The instruments and implants in this system are intended to be used together as described in the Surgical Technique Guide.
- The instruments and implants are designed to be used together and any deviations from the use described in the Surgical Technique may lead to surgical complications and implant failure.

MRI SAFETY INFORMATION

Many Acumed implants have been evaluated for safety in the MR environment. Consult our publication "Acumed Implants in the MR Environment" at www.acumed.net/ifu for more information.

STERILITY

- Implants and instruments may be provided either sterile or non-sterile as indicated on the label.
- Non-sterile devices are intended to be sterilized before use.

IMPLANTS

MATERIALS

- The implant is made of commercially pure Titanium per ASTM F67.
- Consult our document "Metal Sensitivity Statement" at www.acumed.net/ifu for the chemical composition of Acumed metal implants.
- Chemical characterization data confirmed that all identified chemicals and elements were present at levels below 0.1 w/w%. Toxicological assessment further showed that all identified values were below tolerable exposure levels and are biologically safe.

SINGLE USE

- Implants are intended for single use only, as indicated on the label.
- Do not reuse single use implants as this may increase the risks of failure and cross-contamination.
- Dispose of any unused implant that is contaminated with human blood or tissue. Do not process a contaminated implant.

DEVICE HANDLING

- For safe and effective use, the surgeon must be thoroughly familiar with the implant, the methods of application, instruments, and the recommended surgical technique.
- Physiological dimensions limit implant sizes. Select the type and size of implant that best meets the patient's requirements for close adaptation and firm seating with adequate support.
- Implants are designed to hold fractured bones in place during bone healing at the site of fracture. Implants are not designed for excessive loadbearing. Improper selection or improper implantation of the device may increase the possibility of loosening or migration.
- Do not bend the implant except as indicated in the surgical technique. Repeated or excessive bending may weaken the implant and cause failure at a later time. Refer to the surgical technique.
- Only combine implants when they are intended for that purpose.
- Protect implants against scratching and nicking to prevent stress concentrations, which can result in failure.
- Prevent unused implants from becoming soiled.
- The color of anodized implants may change over time due to processing. This color change does not affect the mechanical properties of the implants.

SAFE DISPOSAL

Implant components and single use, non-powered instruments are to be disposed of in accordance with applicable laws, rules, and regulations for the disposal of biohazardous waste. Follow all guidelines for biohazardous waste in accordance with the Centers for Disease Control and Prevention guidelines as well as applicable federal/national, state and local regulations.

Please follow specific country, state/province and local instructions for safe disposal of electrical devices.

INSTRUMENTS

MATERIALS

The instruments are manufactured from various grades of titanium, stainless steel, aluminum and Radel.

MULTIPLE USE and SINGLE USE

- Instruments are intended for multiple use unless identified on the label for single use only.
- Single use instruments are intended to be disposed after use on a single patient during a single procedure.
- Do not reuse single use instruments as this may increase the risks of failure and cross-contamination.
- Multiple use instruments are only intended for use on a single patient and a single procedure before requiring processing.

- Multiple (limited) use instruments, have a limited lifespan. Immediately replace any multiple use instrument if performance becomes inadequate.
- Multiple use instruments potentially contaminated with transmissible spongiform encephalopathy (TSE) agents shall not be processed for reuse.

IMPORTANT

- Protect instruments against scratching and nicking to prevent stress concentrations, which can lead to instrument failure.
- Near the point of use: Wipe excess contamination from instruments and prevent any soil from drying. Instruments with substantial or dried soil are particularly difficult to reliably process. Transport contaminated instruments for processing as soon as possible after use.
- Avoid prolonged instrument contact with iodine and saline.
- Handle and transport soiled instruments in a manner that avoids contamination of any unused implants.

PROCESSING

Important: Processing personnel must be qualified with suitable training and experience. Use proper personal protective equipment (PPE) when working with contaminated devices.

IMPORTANT

- Promptly perform the processing steps to limit microbial growth and maximize the effectiveness of sterilization.
- Inspect implants for contamination by blood or tissue and dispose of them when found. Do not process contaminated implants.
- Prevent instrument corrosion by minimizing contact with solutions containing iodine, chlorine, and saline or other metal salts.
- Prevent damage to the protective anodization layer on aluminum instruments by avoiding contact with solutions < 4 pH and > 9 pH, especially if they contain sodium carbonate or sodium hydroxide.
- Repeated processing of anodized metals may cause colors to fade but this does not affect the function of the device.
- Avoid cleaning agents containing aldehydes since they can denature and coagulate proteins (fixation).
- Enzymatic detergents are well suited for loosening protein-based contamination.
 - Use a neutral pH enzymatic detergent.
 - Use a low foaming solution to allow visibility of the device during cleaning.
- Closely follow the manufacturer's instructions for the safety, storage, mixing, water quality, exposure time, temperature, replacement, and disposal of cleaning agents.
- Devices potentially contaminated with transmissible spongiform encephalopathy (TSE) agents shall not be processed or reused. These processing instructions are not suitable for inactivation of TSE agents. Observe hospital procedures, practice guidelines, and/or government regulations for the proper handling and disposal of devices potentially contaminated with TSE agents.
- Utility water: Refer to AAMI TIR34* when instructed to use utility water. Utility water is typically municipal or tap water but may require additional treatment to be suitable for use.
- Critical water: Refer to AAMI TIR34* when instructed to use critical water. Critical water is highly treated and has very low organic and inorganic content with an endotoxin level under 10
- EU/mL. Suitable water may also be specified in national pharmacopeias, national standards, and hospital protocols.

* Association for the Advancement of Medical Instrumentation (AAMI). Water for the reprocessing of medical devices. AAMI TIR34:2014/(R)2017. Arlington, VA

MANUAL CLEANING

1. Rinse the contaminated instruments under running cold utility water to reduce heavy surface contamination.
2. Dispose of any used instruments intended for single use only.
3. Place the contaminated instruments in enzymatic solution* until completely submerged to minimize the spraying of solution.
4. Actuate all moveable parts to allow detergent to contact all surfaces
5. Soak for a minimum of ten (10) minutes.
6. Scrub the instruments using a soft-bristled brush to remove all visible debris. Do not use stainless steel or other abrasives as these may damage the surface.
 - When possible, scrub the instruments when totally submerged to minimize the spraying of fluid.
7. Some instruments may require special consideration:
 - Clean the instruments with all parts loosened. Clean the instruments disassembled if they are designed to be taken apart.
 - Use a water jet to flood cleaning solution into challenging areas, such as mating surfaces, springs, coils, cannulations, blind holes, flutes, cutting teeth, and flexible parts to flush out any trapped soil.
 - Operate movable parts and rotate (as necessary) while scrubbing to ensure that all crevices are accessible.
 - Carefully clean cannulated parts and challenging areas using an appropriately sized brush.
 - Optionally sonicate for 10 to 15 minutes using a fresh, neutral pH ultrasonic cleaning solution. Follow the ultrasonic cleaner and detergent manufacturer's instructions.
Important! Any previous surface damage may increase due to ultrasonic cleaning.
8. Perform an initial rinse for at least 3 minutes using clean, soft, utility water in the temperature range of 25°C to 35°C (77°F to 95°F) to remove all signs of contamination and cleaning agent.
 - Actuate all movable parts.
 - Flush out cannulations and complex mechanisms.
9. Repeat the previous processing steps if visible residue remains present.
10. Perform a final rinse for at least 1 minute using critical water to displace minerals and other impurities found in utility water. Do not use saline solutions for final rinsing because they may interfere with disinfection and sterilization.
 - Actuate all moving parts
 - Pay particular attention to cannulations and blind holes as well as hinges and joints between mating surfaces.
 - Rinse cannulations at least three times with a syringe (volume 1-50ml).
11. Remove excess moisture from the instruments using a clean, absorbent, non-shedding wipe. 12. Allow the instruments to thoroughly dry. Any moisture may affect sterilization and devices may remain wet after the drying period.

* Manual cleaning was validated using STERIS Prolystica 2X Concentrate Enzymatic Presoak and Cleaner.

AUTOMATED CLEANING

Important: Automated cleaning requires manual cleaning. Always perform the previous manual cleaning steps first. Follow the washer-disinfector manufacturer's instructions explicitly.

Equipment: An automated washer-disinfector that has been installed and qualified to ISO 15883-1 and ISO 15883-2.

Solutions: Prepare solutions per manufacturer's instructions. Use Neutral pH (<8.5) low foaming, enzymatic wash solution (e.g., Enzol®), neutral pH, low foaming wash solution (e.g., Prolystica Neutral 2x), instrument lubricant (e.g., Ultra Clean Surgical Milk) if capable.

NOTE: Explicitly follow washer/disinfector manufacturer's instructions for loading.

Motor Speed: High

| Phase | Recirculation Time | Temperature | Detergent Type and Concentration *or equivalent |
|-------------|--------------------|----------------|---|
| Pre-wash | 2.00 minutes | Cold tap water | N/A |
| Enzyme Wash | 4.00 minutes | Hot tap water | *Enzol® 1 oz/gal |
| Wash | 2.00 minutes | 65.5°C (150°F) | *Prolystica 2x Neutral 1/8 oz/gal |
| Rinse | 15 seconds | Hot tap water | N/A |
| Drying | 6.00 minutes | 98.8°C (210°F) | N/A |

PRE-STERILIZATION INSPECTION

- Visually inspect all devices under normal lighting to ensure that cleaning was effective. Pay close attention to all challenging areas.
 - Re-process an instrument that is not clean.
 - Replace an instrument that cannot be cleaned.
- Inspect the implants and instruments for surface damage, such as nicks, scratches, and cracks. Replace any device that is affected.
- Assess the instruments for proper use. Operate all parts and connecting mechanisms. Give careful attention to drivers, drill bits and reamers, and instruments used for cutting or implant insertion. Critically assess them for wear, sharpness, straightness, and corrosion. Replace any instrument that does not perform as intended.
- Inspect all cutting edges under magnification.
 - Replace an instrument if a cutting edge is dull, chipped, cracked, rolled, or otherwise deformed.
 - Running a cotton cloth over the edge may help detect chipping and cracking.
- Verify the legibility of all markings and reference scales. Replace any device that is unreadable.
- Repair, replace, and/or repeat the cleaning of instruments as needed to ensure proper operation before proceeding with sterilization.
- Lubrication ("instrument milk") may increase the useful life of surgical instruments. Do not use silicone-based lubricants, oil, or grease, as these will interfere with steam sterilization. Only use a water-based lubricant intended for use on surgical instruments and with steam sterilization. Use the lubricant as directed by the manufacturer. Use critical water if dilution is required.
- Fully replenish the system trays and caddies.

STERILIZATION

- Perform sterilization using a dynamic-air-removal (prevacuum) autoclave. Gravity displacement sterilization is not recommended. Immediate use (flash) sterilization is not recommended.
- Ensure the sterilizer's maximum load limit is not exceeded when sterilizing multiple sets or devices.
- Do not stack containers as this might prevent the penetration of steam and inhibit drying.

- Refer to the sterilizer manufacturer's instructions and ensure proper installation, calibration, use, and ongoing maintenance.
 - The sterilized items should be allowed to cool to room temperature before handling. This allows for safe handling and preventing condensation.
 - Follow current industry best practice guidelines such as ANSI/AAMI ST79:2017*.
- * Association for the Advancement of Medical Instrumentation (AAMI). Comprehensive guide to steam sterilization and sterility assurance in health care facilities. AAMI ST79:2017. Arlington, VA.
- The following table shows the minimum parameters validated* to achieve a required Sterility Assurance Level (SAL) of 10⁻⁶ for the system.

Important:

- Sterilization parameters are only valid for devices that have been cleaned per these instructions and are thoroughly dry.

Sterilization parameters are only valid when the devices are properly housed in the Acumed storage case part numbers identified in the table.

Prevacuum Steam Sterilizer Parameters

| System Tray Part Number | System Tray Approximate Dimensions | Pre-Vacuum Autoclave | | |
|-------------------------|--------------------------------------|----------------------|------------|----------|
| | | Temperature | Cycle Time | Dry Time |
| RBL4020 | 26 x 54 x 8 cm (10 x 21 x 5 inch) | 132°C (270°F) | 4 min | 30 min |
| | | 135°C (275°F) | 3 min | 25 min |

POST STERILIZATION INSPECTION

- Do not store or use sterile devices if they are not dry.
 - Moisture supports the survival of microorganisms.
 - Moisture remaining on wrapped or contained products after sterilization could compromise the sterile barrier.
 - Moisture can corrode metal and dull sharp edges.

Inspect the sterile barrier for signs of damage. Do not use the product if the sterile barrier has been compromised.

STORAGE CONDITIONS

- Items should be stored under controlled conditions in a manner that minimizes the potential for contamination per ANSI/AAMI ST79:2017.
- Refer to sterilization wrap or rigid container manufacturer's IFU for limits on sterile product storage time and storage requirements for temperature and humidity.

STORAGE OF PACKAGED NON-STERILE

Final packaged product should be stored at room temperature and protected from direct sunlight, pests, and high humidity.

Symbols Glossary

| Symbol | Description | EN ISO 15223-1 |
|---|--|---------------------|
|  | Consult the electronic instructions for use (eIFU) at www.acumed.net/ifu | 5.4.3 |
|  | Caution | 5.4.4 |
|  | Non-sterile | 5.2.7 |
|  | Catalogue number | 5.1.6 |
|  | Batch code | 5.1.5 |
|  | Authorized representative in the European Community / European Union | 5.1.2 |
|  | Medical device | 5.7.7 |
|  | Manufacturer | 5.1.1 |
|  | Date of manufacture | 5.1.3 |
|  | Do not re-use | 5.4.2 |
|  | MR Conditional | ASTM F2503-20 |
|  | Do not use if package is damaged and consult instructions for use / do not use if the product sterile barrier system or its packaging is compromised | 5.2.8 |
| Rx Only | Caution: U.S. federal law restricts this device to sale by or on the order of a physician. | U.S. 21 CFR 801.109 |
|  | The reticle is a registered trademark of Acumed. It may appear alone or with the Acumed name. | |
|  | CE marking of conformity, Article 17 of EU Directive 93/42/EEC or Article 20 of Regulation (EU) 2017/745. CE marking may be accompanied by the identification number of the notified body responsible for conformity assessment. | |

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The prothesis is intended to be press fit or cemented.

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These materials contain information about products that may or may not be available in any particular country or may be available under different trademarks in different countries. The products may be approved or cleared by governmental regulatory organizations for sale or use with different indications or restrictions in different countries. Products may not be approved for use in all countries. Nothing contained in these materials should be construed as a promotion or solicitation for any product or for the use of any product in a particular way that is not authorized under the laws and regulations of the country where the reader is located. Nothing in these materials should be construed as a representation or warranty as to the efficacy or quality of any product, nor the appropriateness of any product to treat any specific condition. Physicians may direct questions about the availability and use of the products described in these materials to their authorized Acumed distributor. Specific questions patients may have about the use of the products described in these materials or the appropriateness for their own conditions should be directed to their own physician.

PKGI-011480-D

| Effective 07-2024 | © 2024 Acumed® LLC

Acumed® is a registered trademark of Acumed, LLC.

Subsidiaries:

Acumed Ltd

Huebner House
The Fairground
Andover
Hampshire UK SP11 0QN
Tel: +44 1264 774450

Acumed Iberica

C/ Proción 1
Edificio Oficor
28023 Madrid, Spain
Tel: +34 913516357

Acumed Beijing

Room A1206, Horizon International
Tower
No. 6, Zhichun Road
Haidian District
100088 Beijing, China
Tel: +86 10 82001303

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Acumed GmbH