


 acumed[®]
Acu-Sinch[®] Knotless Mini System



Bringing Acu-Sinch Knotless
simplicity to smaller, more
delicate **soft tissue repairs**




Drill. Deploy. Sinch.
It's that simple.

Acu-Sinch® Knotless Mini System

The Acu-Sinch Knotless Mini is designed for anatomically challenging, smaller repairs, maximizing strength while minimizing complexity. With a low-profile, surgeon-friendly design for small joint precision, it supports fast procedures, predictable outcomes, and strong fixation in soft tissue repair.

Drill. Deploy. Sinch. It's that simple.



Handle Trigger

Spring-loaded one-step delivery reduces procedural steps and supports efficient implantation

Central Ejector Shaft

Provides smooth, consistent deployment for predictable, seamless small joint repair

Keeled Titanium Flip-Button

Deploys through a 2.3 mm inserter handle for streamlined, minimally invasive procedures



Drill

An all-suture, knotless design enables fast execution, preserves soft tissue, and eliminates knot-tying complexity for a more seamless surgical procedure.

Deploy

Intuitive instrumentation designed to deliver a consistent, predictable, and reliable surgical procedure.



Low-Profile Round Buttons

Less than 1 mm thick to minimize soft tissue irritation and support controlled tensioning




Slotted Drill Design

Maintains tunnel alignment for precise, reliable placement

Sinch

Maintains strong, uniform soft tissue tension without slipping or knot failure.




2-0 UHMWPE Suture

Engineered for durability to provide consistent, reliable soft tissue fixation

Self-Locking Suture Loop

Knotless design secures soft tissue while maintaining construct integrity

Contact Your Local Sales Rep for More Information



Learn More at
go.acumed.net/TL-25-ASK-Mini

Drill. Deploy. Sinch.
It's that simple.



www.acumed.net [in @acumed](#) [ig @acumedglobal](#)

Acumed Oregon Campus
5885 NE Cornelius Pass Road
Hillsboro, OR 97124
+1.888.627.9957

Acumed Texas Campus
3885 Arapaho Road
Addison, TX 75001
+1.800.456.7779

Acumed Iberica Campus
C. Proción, 1
Edificio Oficor
28023 Madrid, Spain
+34.913.51.63.57

TLB10-02-B | Effective: 2025/09 | © 2025 Acumed® LLC

Acumed® and Acu-Sinch® are registered trademarks of Acumed LLC.